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**LDI MS analysis.** Laser desorption/ionization (LDI) time-of-flight (ToF) mass spectrometry experiments were performed using a Bruker Autoflex Speed reflectron time-of-flight mass spectrometer, equipped with a SmartBeam II laser (352 nm) in 80-2000 m/z range. The laser impulse energy was approximately 60–120  $\mu$ J, the laser repetition rate was 1000 Hz, and the deflection value was set on  $m/z < 80$  Da. The first accelerating voltage was held at 19 kV, and the second ion-source voltage was held at 16.7 kV. The reflector voltages used were 21 kV (first) and 9.55 kV (second). The data was recorded and analyzed using the software provided with the Autoflex instrument (FlexAnalysis version 3.3). Mass calibration (typically cubic calibration based on five to seven points) was performed using internal standards (gold ions and clusters from  $Au^+$  to  $Au_{10}^+$  depending on  $m/z$  range).

#### Wykonanie:

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**LDI MS sample preparation.** Sample solution (ca. 5 mg/ml in  $H_2O$ ) was placed on AuNPET (0.5  $\mu$ L). Sum of ca. 7000 scans was collected for each sample.

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**LDI MS sample preparation.** Sample solution (ca. 5 mg/ml in  $H_2O$ ) was placed on AuNPET (0.5  $\mu$ L) with 0.5  $\mu$ L standard CHCA solution (1:1 water:acetonitrile with 0.2% TFA). Sum of ca. 7000 scans was collected for each sample.

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Z komentarzem [t4]: jw

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J. Sekuła, J. Nizioł, W. Rode and T. Ruman, “Gold nanoparticle-enhanced target (AuNPET) as universal solution for laser desorption/ionization mass spectrometry analysis and imaging of low molecular weight compounds”, *Analytica Chimica Acta* 875, 2015, 61–72.

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